



3-Ply Surgical Mask - Level 1

Specifications

- ✓ Outer-Layer Spunbond NWPP (25g), Mid-Layer Meltblown Fiber (25g), Inner-Layer Spunbond NWPP (25g)
- ✓ Highly Breathable and Lightweight with Moldable Nose Strip
- ✓ Latex-free TPU Ear Loops & Ultrasonic Construction
- ✓ FDA Registered Class 2 Medical Device (FXX)
- ✓ FDA 510(k) Premarket Notification

**FDA
510(k)**
Premarket
Notification

[Click to Verify](#)

LEVEL 1 (ASTM F2100) - *Minimum Guaranteed Performance

BFE (3.0 µm)	≥ 95%*	≥95% BFE Bacterial Filtration Efficiency
PFE (0.1 µm)	≥ 95%*	
FLUID RESISTANCE (mmHg)	120	Click To View Report
DELTA P (mm H ₂ O/cm ²)	<6.0	
FLAME SPREAD	Class 1	

Product Details

Model	KF-B P05
Mask Size	6.69" x 3.74"
Standard Box Size	50 count
Boxes per Case	50
Cases per Pallet	30



U.S. Sales & Service:

(800) 674-2894
sales@rocket.direct

KINGFA | MEDICAL
Authorized Distributor

rocket.direct