

3-Ply Surgical Mask - Level 2



Specifications

- ✓ Outer-Layer Spunbond NWPP (25g), Mid-Layer Meltblown Fiber (25g), Inner-Layer Spunbond NWPP (25g)
- ✓ Highly breathable and lightweight with moldable nose strip
- ✓ Latex-free TPU ear loops & ultrasonic construction
- ✓ Provides fluid protection barrier (120 mm Hg)
- ✓ Level 2 (ASTM F2100-11); Type IIR (EN14683)

≥98%
BFE
Bacterial Filtration
Efficiency

Product Details

| | |
|-------------------|--------------|
| Mask Size | 6.69 x 3.74" |
| Standard Box Size | 50 count |
| Boxes per Case | 50 |
| Cases per Pallet | 30 |
| Pallets per Truck | 28 |



Fully Customized Packaging and Brand Solutions:



Individually Wrapped
(Available with Custom Print)



Custom Quantity
Retail Packs



Custom Color and/or
Custom Printing

U.S. Sales & Service:

(800) 674-2894
sales@rocket.direct

rocket.direct