



KN95 Respirator

Specifications

- ✓ Outer-Layer Spunbond NWPP (50g), Mid-Layer Meltblown Fiber (2 x 25g), Inner-Layer Spunbond NWPP (30g)
- ✓ Latex-free TPU ear loops & ultrasonic construction
- ✓ Moldable nose strip for optimal seal and fit
- ✓ Flow rate 85 L/min
- ✓ KN95 (GB2626-2006); FFP2 (EN149-2001)

≥95%
PFE
0.3 Micron Particulate
Filtration Efficiency

Product Details

Mask Size:	Adult
Standard Box Size	50 count
Boxes per Case	50
Cases per Pallet	24
Pallets per Truck	20



Custom Retail Packaging Available



FDA Authorized - Appendix A Listed

Guangdong Kingfa Sci. & Tech. Co. Ltd.
Model Number KF-A F01

U.S. Sales & Service:

(800) 674-2894
sales@rocket.direct

KINGFA
金发科技

Authorized Distributors

rocket.direct