



3-Ply ASTM Level 2 Mask - Youth Size

Specifications

- ✓ Outer-Layer Spunbond NWPP (25g), Mid-Layer Meltblown Fiber (25g), Inner-Layer Spunbond NWPP (25g)
- ✓ Highly breathable and lightweight with moldable nose strip
- ✓ Latex-free TPU ear loops & ultrasonic construction
- ✓ Provides fluid protection barrier (120 mm Hg)
- ✓ Level 2 (ASTM F2100-11); Type IIR (EN14683)

≥98%
BFE
Bacterial Filtration
Efficiency

Product Details

Mask Size	5.7" x 3.7"
Standard Box Size	50 count
Boxes per Case	50
Cases per Pallet	30
Pallets per Truck	28



Fully Customized Packaging and Brand Solutions:



Individually Wrapped
(Available with Custom Print)



Custom Quantity
Retail Packs



Custom Color and/or
Custom Printing

U.S. Sales & Service:

(800) 674-2894
sales@rocketdistributors.com

Rocket DISTRIBUTORS