

3-Ply Surgical Mask - Level 3



Specifications

- ✓ Outer-Layer Spunbond NWPP (25g), Mid-Layer Meltblown Fiber (25g), Inner-Layer Spunbond NWPP (25g)
- ✓ Highly Breathable and Lightweight with Moldable Nose Strip
- ✓ Latex-free TPU Ear Loops & Ultrasonic Construction
- ✓ FDA Registered Class 2 Medical Device (FXX)
- ✓ FDA 510(k) Premarket Notification (K202139)

**FDA
510(k)**
Premarket
Notification

[Click to Verify](#)

LEVEL 3 (ASTM F2100-19) - * *Minimum Guaranteed Performance*

BFE (3.0 μm)	$\geq 98\%^*$	Nelson Labs® Tested $\geq 99.9\%$ BFE Bacterial Filtration Efficiency
PFE (0.1 μm)	$\geq 98\%^*$	
FLUID RESISTANCE (mmHg)	160	Click To View Report
DELTA P (mm H ₂ O/cm ²)	<6.0	
FLAME SPREAD	Class 1	

Product Details

Model	KF-B P05(L3)
Mask Size	6.69" x 3.74"
Standard Box Size	50 count
Boxes per Case	50
Cases per Pallet	30



U.S. Sales & Service:

(800) 674-2894
sales@rocket.direct

KINGFA | MEDICAL
Authorized Distributor

rocket.direct