


4 mil Chemo-Rated Nitrile Examination Gloves - 510(k)



Product Box
100 gloves




Master Carton
10 boxes / case
1,000 gloves



Specifications

- ✔ 100% Nitrile (non-blend)
- ✔ Medical grade
- ✔ Chemotherapy drug resistant
- ✔ Latex and powder free
- ✔ Textured finger
- ✔ Color: Indigo



Size	Median Length (mm)	Median Width (mm)
S	≥ 240	80 ± 10
M	≥ 240	95 ± 10
L	≥ 240	110 ± 10
XL	≥ 240	≥ 110

Technical Details:

Manufacturer:	Kingfa Sci. & Tech. Co., Ltd.
Model Number:	KG-1802
510(k) Number:	K203593
Certifications:	ASTM D6319 & D6978, EN 374, EN 455
Cuff length:	Standard
Glove length:	240 mil (9.5 inches)
Finger thickness:	4.7 ± 0.5 mil
Palm thickness:	3.2 ± 0.5 mil
Freedom from holes:	1.5 AQL
Breakage force:	≥ 6N
Tensile strength:	≥ 14 Mpa
Elongation:	≥ 500%
Shelf life:	3 years

U.S. Sales & Service:

(800) 674-2894
sales@rocketdistributors.com

KINGFA | MEDICAL
Authorized Distributor

Rocket DISTRIBUTORS